



PEROXigard™

One-Step Disinfectant Cleaner & Deodorizer
Wipes

FOR USE IN LIFE SCIENCES

For Use in Life Science Research Facilities, Vivariums, Laboratories, and Animal Research Facilities.

For disinfection and cleaning of hard non-porous surfaces found in animal premises facilities such as biosafety cabinets, lab benches, animal transportation vehicles, animal housing rooms, animal examination rooms, animal imaging and diagnostics rooms, and general laboratory equipment.

Claims

A Powerful Disinfectant Effective Against:

TUBERCULOCIDAL: 3 MIN.

Mycobacterium terrae

VIRUCIDAL: 3 MIN.

Proven effectiveness against the Poliovirus Type 1, Sabin strain type 1 which allows for Broad-Spectrum Virucide claim against most enveloped and non-enveloped viruses.

HIV-1 Human Immunodeficiency virus (HIV), Strain HTLV-IIIB (HIV-1)

Rhinovirus serotype 14

Canine Parvovirus (CPV), the Cornell strain

Feline Calicivirus, F9 Strain, as a surrogate for Norwalk and Norwalk-like viruses

Murine Norovirus type 1 (MNV-1) strain S99

Human Coronavirus 229E

Human Parainfluenza Virus 3

Influenza Virus PR8 strain

Rotavirus Wa strain

Adenovirus type 5

This product has demonstrated effectiveness against Poliovirus and is expected to inactivate all Influenza A viruses including 2009 (H1N1) pandemic Influenza A virus

BACTERICIDAL: 3 MIN.

Staphylococcus aureus

Pseudomonas aeruginosa

Salmonella enterica

Methicillin-resistant Staphylococcus aureus (MRSA)

Vancomycin-resistant Enterococcus faecalis (VRE)

Escherichia coli O157:H7

Enterobacter aerogenes

FUNGICIDAL: 3 MIN.

Trichophyton mentagrophytes

Candida albicans



Reference Sheet

Manufactured by:
Virox® Technologies Inc.
2770 Coventry Road
Oakville, ON Canada L6H 6R1
1-800-387-7578 | virox.com



READ SAFETY DATA SHEET BEFORE USING PRODUCT
SDS available online at www.virox.com

BROAD-SPECTRUM SANITIZING: 30 SEC.

Staphylococcus aureus

Pseudomonas aeruginosa

Salmonella enterica

Klebsiella pneumoniae

Acinetobacter baumannii

Methicillin-resistant Staphylococcus aureus (MRSA)

Vancomycin-resistant Enterococcus faecalis (VRE)

Escherichia coli O157:H7

Campylobacter jejuni

Disinfectant Use Directions

HEAVILY SOILED SURFACES REQUIRE CLEANING PRIOR TO DISINFECTION

Cleaning, Disinfection and Deodorizing of Laboratory Equipment and Hard Non-Porous Surfaces: Pre-clean heavily soiled areas. Apply solution to surface with disposable wipe until thoroughly wet. Ensure surface remains wet for 3 minutes. Allow to air dry.

Directions for Use in Animal Containment Areas

1. Remove all animals and their feed from premises and enclosures prior to disinfection.
2. Remove heavy soil from the surfaces of facilities and fixtures occupied by animals.
3. Empty or remove all feeding and watering equipment.
4. Thoroughly clean all surfaces with detergent or product and rinse with potable water.
5. Apply product to all surfaces for a period of 3 minutes.
6. Allow all equipment used in the handling of animals to remain wet for 3 minutes. Allow items to dry before re-using.
7. Ventilate closed spaces. Do not house animals or employ equipment until the product has been absorbed, set, or dried.
8. Thoroughly scrub all treated surfaces or objects that come in contact with feed and drinking water with soap or detergent, and rinse with potable water before reuse.

Directions for Use in Food Processing Establishments

Apply Peroxigard One-Step Disinfectant Cleaner & Deodorizer Wipes for use in Life Sciences to surface. Ensure surface remains wet for 3 minutes. For surfaces and/or equipment that may come in contact with food, a rinse with potable water is required.

Broad-Spectrum Sanitizing on Environmental Surfaces (>99.9%)

Apply to surface, allow to remain wet for 30 seconds. Allow to air dry. Food contact surfaces require rinsing with potable water.

To Pull Wipes Through Lid, Simply...

1. Remove lid.
2. Remove liner from the canister.
3. Pull up corner of first wipe from centre of roll and push 1 to 2 cm, ½" to ¾" through "X" in lid.
4. Replace lid.
5. Lift lid flip top and pull up wipes one at a time.

To close between uses to prevent drying, simply keep lid flip top down and ensure wipe corner is inside the raised circle in lid.

Precautionary Statements

No Personal Protective Equipment is required when handling this product. **KEEP OUT OF REACH OF CHILDREN.** Avoid storage at elevated temperatures. Avoid contamination of food in the application and storage of this product. Do not mix with other cleaning or disinfecting products.

FIRST AID: If in contact with eyes, flush immediately and thoroughly with water. Call a physician if irritation persists. If ingested, do not induce vomiting. Call a physician or a poison control center immediately. **TAKE CONTAINER, LABEL OR PRODUCT NAME AND DIN WITH YOU WHEN SEEKING MEDICAL ATTENTION.**

STORAGE: Avoid storage at elevated temperatures. Store at room temperature.

DISPOSAL: For information on the disposal of unused, unwanted product and the clean up of spills, contact the Provincial Regulatory Agency or the Manufacturer. Rinse the container with water. Dispose of container in accordance with municipal, provincial and federal regulations. Do not reuse containers.

Special Instructions for Cleaning and Decontamination Against HIV (Human Immunodeficiency Virus) on objects and surfaces soiled with blood/body fluids. This product is intended for use against HIV only in those settings where the virus would be expected to be encountered, such as settings where contamination by blood or body fluids is likely.

Cleaning and Disinfecting Surfaces of Blood and Body Fluids: Gloves should be worn. Remove excess blood and fluid with absorbent materials. Clean contaminated area: Wipe the surface with disposable Peroxigard One-Step Disinfectant Cleaner & Deodorizer Wipes for use in Life Sciences. Ensure all blood/body fluids are thoroughly cleaned from surfaces/objects before starting disinfection. Disinfect contaminated area: Apply to surface with disposable wipe, allow surface to remain wet for 3 minutes. Wipe surface dry or rinse.

Personal Protection: Disposable gloves, gowns, face masks, or eye covering as appropriate, must be worn during all cleaning of body fluids, blood and decontamination procedures.

Disposal of Infectious Material: Products contaminated with blood or body fluids should be disposed of according to Federal, Provincial, and local regulations for infectious waste disposal.

Germicidal activity of this product was determined in accordance with the Canadian General Standards Board's standard CAN/CGSB-2.161-97