

# VTAC51855G AciSolv, Acid Cage Wash Detergent

A versatile low-foaming detergent specifically designed for the removal of inorganic residues and salts.  
Designed for manual and automated cage washers.

## Precautionary Statements:

- Wash hands and exposed areas thoroughly after handling
- Wear eye protection, protective gloves
- If on skin: Wash with plenty of water
- If in eyes: Rinse cautiously with warm water for several minutes
- Use with adequate ventilation. Avoid breathing vapors, spray or mist
- Store in a cool dry place. Keep partially used containers closed and away from heat and humidity. Do not exceed storage temperature of (50 °C) 120°F. Keep out of reach of children.
- Eye protection: Safety glasses are required
- Skin protection: Gloves are required
- Respiratory protection: Use an approved respirator if exposure limits are exceeded
- Hygiene measures: Handle in accordance with good industrial hygiene and safety practice
- Dispose of in accordance with local, state and federal regulations

## In Case of Fire:

- Extinguish Mode: Use extinguishing measures appropriate to local circumstances and surrounding environment.
- Flash Point: Not applicable
- Auto Ignition Temperature: Not applicable
- Specific Hazards: Not flammable or combustible
- Special Protective Equipment: As in any fire, wear a self-contained breathing apparatus (SCBA), approved according to NFPA / NIOSH or equivalent and full protective gear.

## First Aid:

- Eye Contact: Immediately flush with cool water for at least 20 minutes. Seek medical attention immediately.
- Skin Contact: Immediately wash skin with soap and water for at least 20 minutes while removing contaminated clothing. Seek medical attention if irritation persists.
- Ingestion: Do NOT induce vomiting. Drink cool water. Transport to an emergency medical center. Never give anything by mouth to an unconscious person.
- Inhalation: Move victim to fresh air. Seek medical attention immediately.

## Hazard Statements:


- Harmful if swallowed
- Harmful in contact with skin
- Causes serious eye damage

## Directions for Use:

Use 1/4 to 2 ounce per gallon of water depending on the type of equipment used, water conditions, and type of debris.



## WARNING



## Company Info:

- Company Name –Lighthouse Environmental Infection Prevention
- Address 100 TradeCenter Dr, Suite G-700 Woburn, MA 01801USA
- Emergency Phone #: 800.222.1222



**LIGHTHOUSE**  
Environmental Infection Prevention®