



LIGHTHOUSE
Life Sciences

Custom Stainless Steel Vivarium Equipment

Explore our portfolio of best-in-class stainless steel vivarium equipment, or contact us to place a completely customized order to fit your facility's needs.

What if I need specific vivarium equipment that isn't featured?
Don't fret! We've got you covered. We'll work by your side to make sure that your fully customized, stainless steel vivarium equipment meets all of your facility's needs, whatever they may be.

Bulk Transport Truck

- Overall size
28" deep x 60" wide x 70" high
- Constructed of 1-1/4" square x 16-gauge wall tubular frame
- Enclosed on ends/back with 3/8" diameter type 304 stainless steel rods on 3" vertical centers
- Bottom is solid base with 1/2" diameter holes on 3" x 3" centers
- Heavy-duty stainless steel swivel plate casters are 6" diameter x 2" wide



Mobile PPE Supply Cart

- Overall size:
27-1/2" wide x 19" deep x 47" high
- Our Popular Wall-mounted PPE Station, but Mobile!
- Holds all your PPE supplies in one convenient location
- Made of Durable Lab Grade 304 Stainless Steel with #4 brushed polish finish
- All welds full length and ground/polished #4 finish



Feed Bins Carts

- Overall Size: 18-1/2" deep x 26-1/2" wide x 36" high
- Constructed of 1" square x 16-gauge wall tubular frame and solid 18-gauge formed base shelf, all of type 304 stainless steel with #4 brushed polish finish
- Rated for cage wash
- Base shelf features lip formed up along ends and back and slopes to the front for drainage
- Two front clips for storing lid



Autoclavable Sterilizer Rack

- Overall size:
26" W x 35-1/2" L x 59" H
- Constructed of formed 1-1/4" dia x 18-gauge wall tubular frame, 3/8" dia solid end rods on 4" centers, and 16-gauge solid perforated base shelf
- Base shelf perforated with 1/2" dia holes on 3" x 3" centers
- Features four stainless steel swivel plate casters with 5" dia x 1-1/2" wide high-temp autoclavable wheels



Fully Customizable Racks

- Overall size:
24" deep x 36" wide x 60" high
- Constructed of four formed solid 18-gauge shelves and 1" square x 16-gauge wall corner posts, all of type 304 stainless steel with #4 brushed finish. 16-1/2" clearance between shelves (4-shelf model)
- Can also equip with perforated shelves
- All joints welded full length and polished to #4 finish



1-800-298-5433
Your Environmental
Infection Prevention Experts
www.lighthouselifesciences.com

Revision Date: 05/10/22



LIGHTHOUSE
Life Sciences

Custom Stainless Steel Vivarium Equipment

Explore our portfolio of best-in-class stainless steel vivarium equipment, or contact us to place a completely customized order to fit your facility's needs.

Three Shelf Utility Carts

- Overall size: 27" x 40" x 30" high
- Features three 18-gauge formed solid shelves and four 1" square corner posts, all of type 304 stainless steel with #4 brushed finish
- All corners welded full-length and ground/polished to #4 finish
- All shelves have back-edge and both ends formed up 1", 10" clearance between shelves
- All shelves have edge along the back and sides



Replacement Casters And Wheels

- Stainless steel swivel caster with 6" dia x 2" wide glass-filled nylon high-temp autoclavable wheel
- Grease fittings in swivel and axle
- Top plate with 3" x 3" bolt centers, also available with step-on brake
- Full-swivel
- Designed specifically for your customized vivarium equipment
- Replacement casters and wheels for Bulk Truck Transports, Rubbermaid Utility Carts, and Metro Shelving Units



Ready to get in touch?

We'll work side by side with you to create completely customized stainless steel vivarium equipment to meet your facility's needs, whatever they may be.

Or, choose from our portfolio of already designed, best-in-class stainless steel vivarium equipment.

1-800-298-5433
Your Environmental
Infection Prevention Experts
www.lighthouselifesciences.com

Revision Date: 05/10/22