


AciSolv

Section 1: Identification

Product Name:	AciSolv
Product Code:	VTAC51805G, VTAC51830G
Recommended Use:	Low pH detergent
Manufacturer:	Unimed-Lighthouse Distribution
Telephone:	800-347-9023
Website:	www.unimedcorp.com
Issue Date:	18-September-2020
Date of Last Revision:	18-January-2021
Emergency Contact USA:	Phone: 1-800-222-1222
Emergency Contact CA:	Phone: 1-(888)*-226-8832

*666 on a cellular phone

NFPA



Health Hazard	1
Fire Hazard	0
Reactivity	0
Personal Protection	A

Section 2: Hazard(s) Identification

GHS Pictogram:	
Signal Word:	WARNING
Hazard Statement:	H302 - Harmful if swallowed. H312 - Harmful in contact with skin. H318 - Causes serious eye damage.
Precautionary Statements:	P264 - Wash hands and exposed areas thoroughly after handling. P280 - Wear eye protection, protective gloves. P302+P352 - IF ON SKIN: Wash with plenty of water. P305+P351 - IF IN EYES: Rinse cautiously with water for several minutes.

Section 3: Composition / Information on Hazardous Ingredients

Name	CAS No.	% by Weight
Citric Acid	77-92-9	20.0 - 28.0

Section 4: First Aid Measures

Eye Contact:	Immediately flush with cool water for at least 20 minutes. Seek medical attention immediately.
Skin Contact:	Immediately wash skin with soap and water for at least 20 minutes while removing contaminated clothing. Seek medical attention if irritation persists.
Ingestion:	Do NOT induce vomiting. Drink cool water. Transport to an emergency medical center. Never give anything by mouth to an unconscious person.
Inhalation:	Move victim to fresh air. Seek medical attention immediately.

Section 5: Fire Fighting Measures

Extinguishing Media:	Use extinguishing measures appropriate to local circumstances and surrounding environment.
Flash Point:	Not applicable.
Auto Ignition Temperature:	Not applicable.
Specific hazards:	Not flammable or combustible.
Special protective equipment:	As in any fire, wear a self-contained breathing apparatus (SCBA), approved according to NFPA / NIOSH or equivalent and full protective gear.

Section 6: Accidental Release Measures

Personal Precautions:	Keep personnel from entering spill area.
Environmental Precautions:	Avoid runoff into storm sewers, ditches and waterways.
Methods for Containment:	Contain spills with an inert absorbent material, (sand, earth, vermiculite) and place in container for disposal according to local and national regulations.
Methods for Cleanup:	Flush away traces with water and clean with mop or rags.

Section 7: Handling and Storage

Handling:	Use with adequate ventilation. Avoid breathing vapors, spray, or mist.
Storage:	Store in a cool dry place. Keep partially used containers closed and away from heat and humidity. Do not exceed storage temperature of (50 °C) 120°F. Keep out of reach of children.

Section 8: Exposure Controls / Personal Protection

Engineering Controls:	Good general ventilation should be sufficient.
Eye Protection:	Safety glasses are required.
Skin Protection:	Gloves are required.
Respiratory Protection:	Use an approved respirator if exposure limits are exceeded.
Hygiene Measures:	Handle in accordance with good industrial hygiene and safety practice.

Section 9: Physical and Chemical Properties

Appearance:	Clear Liquid	Boiling Point:	> 212°F
Color:	Colorless-Straw Color	Flash Point:	Not determined.
Odor:	Mild	Water Solubility:	Completely soluble.
pH:	< 1.0	Vapor Density:	Not determined.
Specific Gravity:	1.050 -1.10	Vapor Pressure:	Not determined.
Freezing Point:	Not determined.	Evaporation Rate:	Not determined.

Section 10: Stability and Reactivity

Chemical Stability:	Stable under normal conditions.
Hazardous Reactions:	None known.
Conditions to Avoid:	Extreme heat.
Incompatible Materials:	Strong base, nitrocarbons and chlorocarbons.

Section 11: Toxicological Information

Likely Routes of Exposure:	Inhalation, eye contact and skin contact.
----------------------------	---

Potential Health Effects

Eyes:	Causes serious eye irritation.
Skin:	May cause skin irritation. Prolonged contact may cause dryness of the skin.
Inhalation:	Prolonged inhalation of high concentrations may damage respiratory system.
Ingestion:	Gastrointestinal symptoms, including upset stomach. Nausea, vomiting.
Chronic Exposure:	Health injuries are not known or expected under normal use.

Section 12: Ecological Information

Ecotoxicity:	Not determined.
Biodegradability:	Not determined.
Bioaccumulation:	Not determined.

Section 13: Disposal Considerations

Waste Disposal:	Dispose of in accordance with local, state and federal regulations.
-----------------	---

Section 14: Transport Information

Land Transport (DOT):	Not regulated.
Sea Transport (IMDG/IMO):	Not regulated.
Air Transport (IATA):	Not regulated.

Section 15: Regulatory Information

TSCA Inventory:	All components are listed.
CERCLA:	Does not contain any components with reportable quantities.
SARA 302:	None.
SARA 311/312 Hazards:	Acute health hazard.
California Prop 65:	None.

Section 16: Additional Information

Date of preparation:	18-September-2020
Date of Last Revision:	18-January-2021
SDS Author:	Unimed-Lighthouse Distribution

Information contained herein was obtained from sources considered technically accurate and reliable. While every effort has been made to ensure full disclosure of product accurate and reliable. While every effort has been made to ensure full disclosure of product hazards, in some cases data is not available and is so stated. Since conditions of actual product use are beyond control of the supplier, it is assumed that users of this material have been fully trained according to the requirement of all applicable legislation and regulatory instruments. No warranty, expressed or implied, is made and manufacturer/supplier will not be liable for any losses, injuries or consequential damages which may result from the use of or reliance on any information contained in this document. The contents of this document have been prepared in accordance with the OSHA Hazard Communication Standards (2012) and GHS (Globally Harmonized System of Classification and Labeling of Chemicals)