

Who's Guarding You & Your Research?

Peroxigard™ Concentrate for use in Life Sciences

One technology for all your facility's cleaning and disinfection needs. Remains effective for 30 days once diluted. No activation required.

Efficacy

Cleans and achieves disinfection contact time compliance against pathogens of concern in one step. Two times faster than leading concentrate disinfectant products.

Safety

Designed with animal safety in mind. Fragrance free, non-irritating, non-corrosive and superior metal compatibility at the in-use dilution.

Simplicity

Eliminates multi-product programs and activation time, improving productivity and compliance.



DISINFECTION COMPLIANCE
FOR BETTER SCIENCE
Peroxigard.ca



One-Step Disinfectant Cleaner & Deodorizer formulated with Accelerated Hydrogen Peroxide®

Peroxigard™ Concentrate utilizes Accelerated Hydrogen Peroxide®, a patented and synergistic blend of commonly used ingredients, that when combined with low levels of hydrogen peroxide dramatically increase its efficacy and cleaning performance.

Peroxigard™ Concentrate is a simple and effective way to clean and disinfect floors, ceilings and walls within your life science research setting. Our revolutionary technology streamlines your process with no multipart mixing, no activation and has a long shelf life once diluted reducing waste.

Effective against the following pathogens* in as little as 5 minutes



Bactericidal

Klebsiella pneumoniae; Methicillin-resistant Staphylococcus aureus (MRSA); Salmonella enterica; Pseudomonas aeruginosa

Virucidal

Murine Norovirus; HIV Type 1; Rotavirus; Poliovirus Type 1

Tuberculocidal

Mycobacterium terrae

Fungicidal

Candida albicans; Trichophyton mentagophytes

FOR USE IN:

- Vivariums
- Breeding facilities
- Lab Animal Research facilities
- Life Science laboratories
- Animal Examination rooms

FOR USE ON:

- Floors, ceilings, and walls
- Lab benches and exam tables
- Cages, crates, and carts
- Animal housing areas

SHELF LIFE

- 2 years (Opened/Unopened)
Remains effective for up to 30 days once diluted

USE DILUTIONS

- Cleaning & Sanitizing: 1:128
- Disinfection : 1:40

AVAILABLE FORMATS

- 4L, 20L, 200L (DIN 02470918)

* See reference sheet for a complete list of efficacy claims, additional features and directions for use



Peroxigard™ products are formulated with the patented technology known as Accelerated Hydrogen Peroxide®/AHP®. The benefits and efficacy of AHP® are Health Canada Registered and have been validated by more than 55 peer reviewed clinical studies conducted by scientific organizations and recognized third party researchers.

To learn more about Accelerated Hydrogen Peroxide® or download copies of the peer reviewed studies, please visit: Peroxigard.ca

**DISINFECTION COMPLIANCE
FOR BETTER SCIENCE**

