

Who's Guarding You & Your Research?

Peroxigard™ Ready to Use Liquid and Wipes for use in Life Sciences

One technology for all your facility's cleaning and disinfection needs. No dilution or activation required.

Efficacy

Cleans and achieves disinfection contact time in 1 minute* against pathogens of concern in a single application.

Safety

Eliminates the need to interact with hazardous concentrates. Designed with animal safety in mind. Fragrance free, non-irritating, non-corrosive and superior metal compatibility.

Simplicity

Removes the risks associated with diluting concentrates and eliminates the complexity of using multiple products at various dilution ratios.



DISINFECTION COMPLIANCE
FOR BETTER SCIENCE
Peroxigard.com



One-Step Disinfectant Cleaner & Deodorizer formulated with Accelerated Hydrogen Peroxide®

Peroxigard™ Ready to Use Liquid and Wipes utilize Accelerated Hydrogen Peroxide®, a patented and synergistic blend of commonly used ingredients, that when combined with low levels of hydrogen peroxide create an ideal one-step cleaner and disinfectant.

Peroxigard™ RTU and Wipes disinfect in only 1 minute*, removing the risk of infection associated with improperly diluting and not achieving contact time compliance. The wipe product is a convenient way to remove variables by providing a consistent amount of liquid in every wipe.

Effective against the following pathogens* in 1 minute or less



Bactericidal

Klebsiella pneumoniae; Bordetella bronchiseptica; Methicillin-resistant Staphylococcus aureus (MRSA); Salmonella enterica

Virucidal

Murine Parvovirus; Murine Norovirus; HIV Type 1; Rotavirus; Murine Hepatitis virus; Adenovirus

Tuberculocidal

Mycobacterium bovis

Fungicidal

Candida albicans; Microsporum canis; Trichophyton interdigitale

Peroxigard™ Ready to Use Liquid and Wipes are proven to remain effective in the presence of a 5% soil load, providing disinfection confidence.

FOR USE IN:

- Vivariums
- Breeding facilities
- Lab animal research facilities
- Laboratories

FOR USE ON:

- Biosafety cabinets
- Lab benches and exam tables
- Diagnostic and laboratory equipment
- Cages, crates, and carts

SHELF LIFE

- 2 years (Opened/Unopened)

AVAILABLE FORMATS

- Ready-To-Use Liquid (EPA Reg No. 74559-9)
32 oz, 1 gal, 5 gal, 30 gal, 55 gal
- Presaturated Wipes (EPA Reg No. 74559-10)
160 ct. 6x7

To learn more about Accelerated Hydrogen Peroxide® or download copies of the peer reviewed studies, please visit: Peroxigard.com

* See reference sheet for a complete list of efficacy claims, additional features and directions for use



Peroxigard™ products are formulated with the patented technology known as Accelerated Hydrogen Peroxide®/AHP®. The benefits and efficacy of AHP® are EPA approved and have been validated by more than 45 peer reviewed clinical studies conducted by scientific organizations and recognized third party researchers.

DISINFECTION COMPLIANCE
FOR BETTER SCIENCE

