



**LIGHTHOUSE**  
Life Sciences

**ACI Solv**

**ACID DETERGENT CAGE WASH**

For use in both automated cage washers, manual use on a variety of equipment and ultrasonic units within the life science, veterinary, and research industries.



- Versatile, low-foaming detergent specifically designed for the removal of inorganic residues and salts
- Effectively removes mineral scale, hard water deposits, animal residue, and dried soil
- Enhances cleaning and reduces rusting and corrosion
- Can be used in conjunction with **Alka Gent** or other alkaline detergents to equalize residual pH
- Buffer system provides two modes of cleaning power
- Reduces bio-burden from depositing on equipment

## ORDERING INFO

ITEM NUMBER	DESCRIPTION	FORMAT
VTAC51805G	ACI Solv Acid Detergent Cage Wash	5 Gallon
VTAC51830G	ACI Solv Acid Detergent Cage Wash	30 Gallon
VTAC51855G	ACI Solv Acid Detergent Cage Wash	55 Gallon